

**Determination of pyrrolizidine alkaloids (PA)
in honey by SPE-LC-MS/MS**

Method Protocol

BfR-PA-Honey-1.0/2013

CONTENT

1	Scope	1
2	Principles	1
3	Reagents	1
3.1	General	1
3.2	Chemicals	2
3.3	Solutions	2
3.3.1	Extraction solution: 0.05 M sulphuric acid.....	2
3.3.2	Elution solution for solid phase extraction (SPE)	2
3.3.3	HPLC mobile phase	3
3.3.4	Standard solution for calibration	3
4	Apparatus.....	4
5	Procedure.....	5
5.1	Extraction (Betteridge et al. 2005; modified).....	5
5.2	SPE procedure.....	5
5.3	Reconstitution of the sample	5
6	HPLC-MS/MS.....	5
6.1	Liquid chromatographic separation	5
6.2	Mass spectrometric operating conditions	6
6.3	Measurement.....	6
7	Calculation	6
7.1	Calibration function	6
7.2	Quantification	7
7.3	Reporting results.....	7
8	Annex.....	8
8.1	LC-MS/MS measurement.....	8
8.2	Typical chromatogram.....	11
8.3	Provider of PA standards	12
8.4	Flow chart of the sample preparation procedure	15

1 Scope

Pyrrolizidine alkaloids (PA) are secondary plant metabolites with carcinogenic and genotoxic properties. Currently, more than 600 PA are known. They occur in plants of the families of Boraginaceae, Asteraceae and Fabaceae. The pollen of these plants are a potential source for pyrrolizidine alkaloids and may lead to a contamination of honey (EFSA Panel on Contaminants in the Food Chain (CONTAM) 2011).

This method describes the determination of the following PA in honey: senecionine (Sc), senecionine-N-oxide (ScN), seneciphylline (Sp), seneciphylline-N-oxide (SpN), monocrotaline (Mc), monocrotaline-N-oxide (McN), retrorsine (Re), heliotrine (Hn), heliotrine-N-oxide (HnN), trichodesmine (Td), retrorsine-N-oxide (ReN), echimidine (Em), intermedine (Im), lycopsamine (La), senkirkine (Sk), lasiocarpine (Lc) and lasiocarpine-N-oxide (LcN).

The limits of quantification for toxins determined during in-house validation are listed in the Annex (8.1).

2 Principles

A test portion is dissolved in aqueous sulphuric acid solution for PA extraction. Dissolved samples are centrifuged to remove particles and afterwards purified by a solid-phase extraction (SPE) using strong cation exchange material. PAs are released from the cartridge using ammoniacal methanol. Subsequently, the eluate is evaporated to dryness and reconstituted in methanol/water (initial HPLC conditions).

For chromatographic separation, an RP-HPLC column is used with a binary gradient. Analytes are detected by triple stage quadrupole mass spectrometry. Quantification of pyrrolizidine alkaloids is accomplished by means of matrix matched calibration.

3 Reagents

3.1 General

Please note: Since the use of this method involves reagents harmful to health, appropriate precautionary and protective measures such as avoiding skin contact and using an extractor hood must be taken.

If not specified otherwise, reagents of analytical grade and solvents suitable for HPLC-MS/MS must be used. Water must be distilled in glass vessels or demineralised before use, or must be of equivalent purity.

3.2 Chemicals

- 3.2.1 echimidine (Em)
- 3.2.2 heliotrine (Hn)
- 3.2.3 heliotrine-N-oxide (HnN)
- 3.2.4 intermedine (Im)
- 3.2.5 lasiocarpine (Lc)
- 3.2.6 lasiocarpine-N-oxide (LcN)
- 3.2.7 lycopsamine (La)
- 3.2.8 monocrotaline (Mc)
- 3.2.9 monocrotaline-N-oxide (McN)
- 3.2.10 retrorsine (Re)
- 3.2.11 retrorsine-N-oxide (ReN)
- 3.2.12 senecionine (Sc)
- 3.2.13 senecionine-N-oxide (ScN)
- 3.2.14 seneciphylline (Sp)
- 3.2.15 seneciphylline-N-oxide (SpN)
- 3.2.16 senkirkine (Sk)
- 3.2.17 trichodesmine (Td)
- 3.2.18 formic acid 98 – 100%, e.g. Sigma-Aldrich
- 3.2.19 methanol (MeOH) in LC-MS quality, e.g. Merck LiChrosolv®
- 3.2.20 sulphuric acid 98%, e.g. Merck
- 3.2.21 ammonia 32%, e.g. Merck
- 3.2.22 ammonium formate in LC-MS quality, e.g. Fluka
- 3.2.23 acetonitrile, e.g. Merck LiChrosolv®

3.3 Solutions

- 3.3.1 Extraction solution: 0.05 M sulphuric acid
2.665 mL of sulphuric acid (H₂SO₄) (3.2.20) are filled up to 1 L with water. The final concentration is 0.05 M.
- 3.3.2 Elution solution for solid phase extraction (SPE)
To prepare a 2.5 % ammoniacal solution for the elution of PA, 7.8 mL ammonia (3.2.21) are filled up to 100 mL with methanol (3.2.19).

Please note: The solution must be prepared every working day.

3.3.3 HPLC mobile phase

Eluent A:

315 mg ammonium formate (3.2.22) are dissolved in 5 mL of water, 1 mL of formic acid (3.2.18) is added and filled up to 1 L with water.

Eluent B

315 mg ammonium formate (3.2.22) are dissolved in 5 mL of water, 1 mL of formic acid (3.2.18) is added and filled up to 1 L with methanol (3.2.19).

3.3.4 Standard solution for calibration

Stock solution (0.1 mg/mL):

To create a stock solution, 1 mg of a pyrrolizidine alkaloid standards are weighed in a volumetric flask using an analytical balance (4.3) and filled up to 10 mL with acetonitrile (3.2.23). The concentration of the stock solution is 0.1 mg/mL.

Standard working solution (PA mixture, 1 µg/mL)

For preparation of the standard working solution, respective volumes of each PA stock solution (0.1 mg/mL) are pipetted into a volumetric flask and filled up with acetonitrile (3.2.23), to achieve a concentration of 1 µg/mL.

Preparation of matrix matched standards (MMS)

For a correction of matrix effects a matrix matched calibration is used. In order to obtain the same matrix strength for MMS as samples the blank honey has to be processed as described in section 5. Afterwards, prepare MMS according to Table 1.

Table 1: Matrix matched standards

	Final PA mass concentration in calibration solution	Final PA mass concentration	Aliquot taken from	Aliquot Volume	Aliquot taken from blank honey extract
	ng/mL	µg/kg		µL	µL
MMS_1	1.0	0.1	MMS_3	20	180
MMS_2	5.0	0.5	MMS_5	20	180
MMS_3	10.0	1.0	MMS_5	40	160
MMS_4	25.0	2.5	Standard working solution (1 µg/mL)	5	195
MMS_5	50.0	5.0	Standard working solution (1 µg/mL)	10	190
MMS_6	75.0	7.5	Standard working solution (1 µg/mL)	15	185
MMS_7	100.0	10.0	Standard working solution (1 µg/mL)	20	180

4 Apparatus

4.1 General

Usual laboratory glassware and equipment should be used and, in particular, the following:

- 4.2 various piston pipettes and multiple dispensers**, e.g. Brand
- 4.3 analytical balance**, capable of weighing to 0,0001 g
- 4.4 compartment drier**
- 4.5 centrifuge** for 50 mL centrifuge tubes, 5000 x g
- 4.6 laboratory shaker**, e.g. Vortex
- 4.7 overhead shaker**, e.g. Heidolph
- 4.8 evaporation station**, e.g. TurboVap
- 4.9 centrifuge tube** 50 mL
- 4.10 test tubes** 15 mL
- 4.11 volumetric flasks**, 10 and 20 mL
- 4.12 SPE cartridges**: DSC-C18 SPE (Supelco), 500 mg sorbent material, 6mL
- 4.13 SPE vacuum chamber**
- 4.14 Membrane filter** 0.2 μm , e.g. VWR 0,5 mL centrifugal filters, modified nylon membrane
- 4.15 HPLC vials** 2 mL
- 4.16 Glass inserts**, 250 μL conic for HPLC vials
- 4.17 chromatographic column**, e.g. Thermo, Hypersil Gold®; 150 x 2.1 mm; 1,9 μm
- 4.18 LC-MS/MS system**

5 Procedure

5.1 Extraction (Betteridge et al. 2005; modified)

For the extraction of PA 10.0 g \pm 0.1 g of honey are weighed into a centrifuge tube (4.9).

The sample is dissolved in 30 mL of extraction solution (3.3.1) at room temperature by shaking (4.7) for 30 min. Sample solutions are centrifuged at 3,800 x g for 10 min \pm 2 min (4.5). The supernatant is transferred to a clean vessel and afterwards used for SPE.

5.2 SPE procedure

SPE (4.12) is carried out using a vacuum chamber (4.13) which is placed in a compartment drier (4.4) at 40 °C \pm 5 °C. By means of the elevated temperature the crystallization of honey at SPE cartridges can be avoided.

Conditioning step 1	5 mL of methanol (3.2.19)
Conditioning step 2	5 mL 0.05 M H ₂ SO ₄ (3.3.1)
Sample load	30 mL sample (supernatant), rinse vessel with 5 mL 0.05 M H ₂ SO ₄ (3.3.1) and add solution to the cartridge
Washing step 1	6 mL of water
Washing step 2	6 mL of methanol (3.2.19)
Drying of cartridges	5 - 10 min (use the vacuum chamber (4.13))
Elution of PA	2 x 5 mL 2.5% ammonia in methanol (3.3.2)

The eluate is dried under a nitrogen stream at 50 °C \pm 5 °C.

5.3 Reconstitution of the sample

The residue is dissolved in 1 mL of methanol/water (5/95, v/v) by shaking (4.6).

The reconstituted sample extracts are filtered through a 0.2 μ m membrane filter (4.14). When using centrifugal filters, 500 μ L of the sample are centrifuged at 20,000 x g for 10 min \pm 3 min. 200 μ L of the filtrate are transferred into an HPLC vial (4.15) with a glass insert (4.16).

6 HPLC-MS/MS

6.1 Liquid chromatographic separation

The measurements can be carried out with different high-performance liquid chromatographs (HPLC) and separation columns. The chromatographic conditions can be chosen freely. The acceptable minimum retention time is twice the retention time for the dead volume of the column. Analytes which cannot be distinguished by means of mass spectrometry must be separated chromatographically. The conditions listed in the annex (8.1) using a C-

18 column (4.17) and the mobile phase described in 3.3.3 have shown to be suitable in pre-trials. However, they are to be seen as examples only.

Please note: For the qualitative detection and for quantification, it is necessary to detect and report two substance-specific transitions per analyte.

6.2 Mass spectrometric operating conditions

The measurements can be carried out with MS/MS devices of different manufacturers. In the annex 8.1, the device-specific settings of one measurement system are given as an example. These conditions have shown to be suitable in pre-trials.

6.3 Measurement

For a quantitative analysis, the following criteria are defined.

Injection:

Samples and standards are injected in duplicate in order to assess repeatability of MS detection and to check for possible response drift during the sequence.

Sequence

To determine pyrrolizidine alkaloids, the following array of analysis is defined in a sequence.

1. Matrix matched standards (1 – 100 ng/mL)
2. Solvent blank
3. Samples (first injection)
4. Solvent blank
5. Matrix matched standards (1 – 100 ng/mL)
6. Solvent blank
7. Samples (second injection)

7 Calculation

The quantitative determination is performed according to the method of the matrix matched standard by integration of the peak areas in relation to the calibration line.

7.1 Calibration function

Equation 1: Calibration function

$$f_{(x)} = y = ax + b$$

where

- | | |
|-----|---|
| y | is the peak area of the target analyte |
| a | is the slope of the calibration function |
| x | is the concentration of the target analyte [ng/mL] in the MMS |
| b | is the intercept of the calibration function |

7.2 Quantification

Equation 2: Calculation of the PA content (analysis equation)

$$PA \text{ concentration} = \beta \times DF = \left[(y - b) \times \frac{1}{a} \right] \times \frac{V_{\text{Extract}}}{m_{\text{weight}}} \times \frac{1}{V_{\text{Application}}} \times V_{\text{sample}}$$

where

β	is the analyte concentration [ng/mL] in the sample extract
DF	is the conversion factor from ng/mL to $\mu\text{g/kg}$
y	is the peak area of the target analyte
a	is the axis intercept from the matrix calibration
a	is the increase from the matrix calibration
V_{Extract}	is the volume of extraction agent [mL]
$m_{\text{Weight-in quantity}}$	is the sample weight in [g]
$V_{\text{Application}}$	is the volume of the extract applied for SPE [mL]
V_{Sample}	is the final sample volume [mL]

7.3 Reporting results

The results are reported in $\mu\text{g/kg}$ with two significant decimals. To convert the concentration from ng/mL injected solution to $\mu\text{g/kg}$ honey a factor of 0.1 is used according to the sample preparation procedure described in chapter 5.

Reference list

- Betteridge K, Cao Y, Colegate SM (2005) Improved method for extraction and LC-MS analysis of pyrrolizidine alkaloids and their N-oxides in honey: application to Echinium vulgare honeys. J Agric Food Chem 53 (6):1894-1902.
- DIN ISO 32645. (1994) Chemical Analysis; Decision limit, Detection limit and determination limit, Estimation in case of repeatability, terms, methods, evaluation. Deutsches Institut für Normung DIN.
- EFSA Panel on Contaminants in the Food Chain (CONTAM). (2011) Scientific Opinion on Pyrrolizidine alkaloids in food and feed. The EFSA Journal 9, 1-135

8 Annex

8.1 LC-MS/MS measurement

LC-MS/MS system consisting of

Triple quadrupole mass spectrometer (Thermo TSQ Vantage)

HPLC system HPLC pump (Thermo Accela 1250),
Degasser
Autosampler (CTC Analytics PAL ATS MYX)
Column oven (MayLab MistraSwitch)

HPLC settings

Eluent A refer to (3.3.3)
Eluent B refer to (3.3.3)
Column temperature 40 °C
Flow rate 300 µL/min
Injection volume 10 µL
Column Thermo Hypersil Gold; 150 x 2.1 mm, 1.9 µm
Total runtime 15 minutes

Gradient

Time (min)	% A	% B
0.0	95	5
0.5	95	5
7.0	50	50
7.5	20	80
7.6	0	100
9.0	0	100
9.1	95	5
15.0	95	5

MS settings

Ionisation Electrospray positive (ESI +)
Ion spray voltage [V] 3500 (positive polarity)
Capillary temperature [°C] 270
Vaporiser temperature [°C] 300
Sheath gas pressure [psi] 45.0
Ion sweep gas pressure [psi] 2.0
Aux gas pressure [psi] 10

Substance-specific parameters

The analytes are detected by Selected Reaction Monitoring (SRM). For analyte identification, two PA specific transitions to two product ions are chosen. The relevant transitions and the collision energy (CE) can be found in Table 2. The table also lists the retention time per analyte which apply for the HPLC settings described above.

Table 2: Substance-specific parameters of the LC-MS/MS method

Analyte	Precursor	Fragment	CE	S Lense	Retention time [min]
Monocrotaline	326.2	120.3	35	130	4,41
		237.3	25	130	
Monocrotaline-N-oxide	342.1	118.3	37	141	5,18
		137.4	29	141	
Intermedine	300.1	138.3	18	112	5,66
		156.3	28	112	
Lycopsamine	300.1	138.3	18	112	5,78
		156.3	28	112	
Retrorsine	352.2	120.3	27	140	6,60
		138.3	29	140	
Trichodesmine	354.2	120.3	35	137	6,68
		222.3	28	137	
Retrorsine-N-oxide	368.2	136.2	30	145	6,71
		118.2	40	145	
Seneciphylline	334.2	120.3	26	138	6,89
		138.4	28	138	
Heliotrine	314.2	138.3	19	119	7,08
		156.3	28	119	
Seneciphylline-N-oxide	350.2	118.2	36	135	7,12
		136.3	32	135	
Heliotrine-N-oxide	330.2	138.2	22	121	7,42
		172.1	27	121	
Senecionine	336.2	120.2	27	130	7,77
		138.2	29	130	
Senecionine-N-oxide	352.2	118.1	28	116	7,98
		136.3	27	116	
Echimidine	398.2	120.3	23	139	8,53
		220.3	17	139	
Senkirkine	366.2	150.3	24	132	8,71
		168.2	28	132	
Lasiocarpine	412.2	120.2	30	139	9,50
		336.3	17	139	
Lasiocarpine-N-oxide	428.2	136.1	29	135	9,75
		254.1	27	135	

Table 3: Limits of detection (LOD) and limits of quantification (LOQ) determined during in-house validation of the described method*

Analyte	LOD [$\mu\text{g}/\text{kg}$]	LOQ [$\mu\text{g}/\text{kg}$]
Monocrotaline	0.22	0.59
Monocrotaline-N-oxide	0.20	0.55
Intermedine	0.13	0.35
Lycopsamine	0.09	0.24
Retrorsine	0.13	0.35
Trichodesmine	0.23	0.62
Retrorsine-N-oxide	0.11	0.29
Seneciophylline	0.10	0.27
Heliotrine	0.13	0.37
Seneciophylline-N-oxide	0.17	0.47
Heliotrine-N-oxide	0.11	0.28
Senecionine	0.15	0.42
Senecionine-N-oxide	0.13	0.35
Echimidine	0.15	0.42
Senkirkine	0.10	0.27
Lasiocarpine	0.07	0.21
Lasiocarpine-N-oxide	0.06	0.18

* LOD and LOQ were determined according to DIN EN ISO 32645 Calibration method (DIN ISO 32645 1994)

8.2 Typical chromatogram




figure 1: Typical chromatogram of a PA standard mixture (10 ng/mL), TIC of SRM-transitions

8.3 Provider of PA standards


pyrrolizidine alkaloid	mass	CAS	provider	order code
echimidine	397,47	520-68-3	Oskar Tropitzsch	7550006
			PhytoLab*	89553
			PlantaAnalytica	-
erucifoline	349.38	40158-95-0	Carl Roth	1657.1
			Oskar Tropitzsch	7550021
erucifoline-N-oxide	365,37	123864-94-8	PhytoLab*	83446
			Carl Roth	1664.1
europine-hydrochloride	365,86	570-19-4	PhytoLab*	83434
			Carl Roth	1676.1
europine-N-oxide	345,39	65582-53-8	PhytoLab*	83237
			AppliChem	A9574,0010
			Carl Roth	1687.1
			Oskar Tropitzsch	7500063
heliotrine	313,40	303-33-3	PhytoLab*	83238
			AppliChem	A9583,0020
			Carl Roth	1929.1
			Latoxan*	L6007
heliotrin-N-oxide	329,39	6209-65-0	Oskar Tropitzsch	7550511
			PhytoLab	80403
			AppliChem	A9590,0010
			Carl Roth	1944.1
indicine-hydrochloride	335,83	1195140-94-3	Oskar Tropitzsch*	755054
			PhytoLab	83236
			Carl Roth	1960.1
indicine-N-oxide	315,36	41708-76-3	Oskar Tropitzsch	7500069
			PhytoLab	83234
			AppliChem	A9593,0010
			Carl Roth	1961.1
intermedine	299,37	10285-06-0	Oskar Tropitzsch	7500070
			PhytoLab	83235
			Carl Roth	1962.1
intermedine-N-oxide	315,36	95462-14-9	Oskar Tropitzsch	7501610
			PhytoLab*	82424
lasiocarpine	411,49	303-34-4	PhytoLab*	83434
			AppliChem	A9596,0010
			Carl Roth	2090.1
			Oskar Tropitzsch*	7500019
lasiocarpine-N-oxide	457,5	127-30-0	PhytoLab	80412
			AppliChem	A9600,0010
			Carl Roth	2202.1
			Oskar Tropitzsch*	7501284
lycopsamine	299,37	10285-07-1	PhytoLab	83220
			Carl Roth	2208.1
			Oskar Tropitzsch	7501080
			PhytoLab*	89726

pyrrolizidine alkaloid	mass	CAS	provider	order code
lycopsamine-N-oxide	315,36	95462-15-0	Oskar Tropitzsch	7501358
			PhytoLab*	83447
monocrotaline	325,35	315-22-0	Carl Roth	3418.1
			Fluka	37024
			Sigma	C2401
			Oskar Tropitzsch	7550522
			PhytoLab*	89251
			R&D Chemicals	7351
monocrotaline-N-oxide	341,36	35337-98-5	Santa Cruz Biotechnology	sc-211921
			Carl Roth	2249.1
			Oskar Tropitzsch	7501658
retrorsine	351,40	480-54-6	PhytoLab*	82629
			AppliChem	A4922,0020
			Carl Roth	1213.1
			Fluka	37025
			Oskar Tropitzsch	7550659
			PhytoLab	89775
retrorsine-N-oxide	367,40	15503-86-3	Santa Cruz Biotechnology	sc-215805
			Sigma*	R0382
			AppliChem	A8668,0010
			Carl Roth	6733.1
senecionine	335,40	130-01-8	Oskar Tropitzsch	7500347
			PhytoLab*	82630
			AppliChem	A2071,0020
			Carl Roth*	2261.1
			Fluka	37031
			Oskar Tropitzsch	7550292
senecionine-N-oxide	351,40	13268-67-2	PhytoLab*	89789
			R&D Chemicals	1828
			Sigma	17806
			Santa Cruz Biotechnology	sc-286770
			AppliChem	A8678,0010
seneciphylline	333,39	480-81-9	Carl Roth	6734.1
			Oskar Tropitzsch	7500301
			PhytoLab*	82631
			AppliChem	A2072,0020
			Carl Roth*	6414.1
			Fluka	37033
seneciphylline-N-oxide	349,38	38710-26-8	R&D Chemicals	1850
			Santa Cruz Biotechnology	sc-229697
			Inc.	
			ABCR GmbH	AB167974
			PhytoLab*	89275
			AppliChem	A8684,0010
			Carl Roth	6735.1
			Oskar Tropitzsch	7500573
			PhytoLab*	82632

pyrrolizidine alkaloid	mass	CAS	provider	order code
senecivernine	335.40	72755-25-0	Carl Roth	2209.1
			Oskar Tropitzsch	7550066
			PhytoLab*	83436
senecivernine-N-oxide	351,39	101687-28-9	Carl Roth	2215.1
			PhytoLab*	83437
senkirkine	365,43	2318-18-5	AppliChem	A6765,0010
			Carl Roth	4934.1
			Fluka	37032
			Oskar Tropitzsch	7500441
			PhytoLab*	89274
trichodesmine	353,41	548-90-3	Latoxan*	L6049

* substances were used for the in-house validation by BfR

8.4 Flow chart of the sample preparation procedure



Bundesinstitut für Risikobewertung